










【Hi-SELECT Containers】

CONTACT

Container dimensions

Product number	Capacity (mL)	Diameter (mm)	Body diameter (mm)	Set height (mm)	Inside diameter (mm)	Materials
No.180	About 180	66.2	61.8	85.2	57.7	Body : HDPE Cap : PP Inner lid : LLDPE
No.300	About 310	83.2	78.7	90.8	74.3	Body : HDPE Cap : PP Inner lid : LLDPE
No.350	About 380	83.2	78.7	105.0	74.3	Body : HDPE Cap : PP Inner lid : LLDPE
No.750	About 800	103.0	98.5	132.3	92.8	Body : HDPE Cap : PP Inner lid : LLDPE
No.1000	About 1,030	103.0	98.5	168.5	92.5	Body : HDPE Cap : PP Inner lid : LLDPE
No.1200	About 1,240	126.6	122.2	132.4	116.2	Body : HDPE Cap : PP Inner lid : LLDPE

Combination

Body	Cap	Inner lid	Notes	Product image
Quantity	Quantity	Quantity		
BHS-180	CSK-150P	SHS-150	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	 
240	360	360		
BHS-300	CSK-350P	SHS-350	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	 
150	300	300		
BHS-350	CSK-350P	SHS-350	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	
150	300	300		
BHS-750	CSK-750P	SHS-750	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	 
80	160	160		
BHS-1000	CSK-750P	SHS-750	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	
40	160	160		
BHS-1200	CSK-1200P	SHS-1200	Cap with inner lid (The inner lid will be set when the cap is snapped on the container.)	
60	120	120		

- The dimensions of No.180 are provisional values.
- Please note that some items may run out of stock.
- Standard colors are the black. Natural color is also available for No.180 and No.300 and No.750 if requested.
- Please note that when the mold is renewed or updated, dimensions and other specifications may change to some degree.
- Air tightness of this product, you may vary depending on content, please review prior to use in the sample. For more information, please contact us.

【Hi-SELECT Containers】



Break the virgin seal in the direction of the mark for one-touch opening.(Tamper-proof mechanism)